



# VitalJacket® SDK


VitalJacket SDK v1.0.07 – BioLib.dll

## LEGAL NOTICE AND DISCLAIMER


**ATTENTION:** Although **VitalJacket** is a certified medical device, its developer version is **NOT** certified for diagnosis usage. It is intended for R&D and development purposes **only**. Users of VJ SDK can submit their final developments to medical certification. All contents of our product are compliant with the European Medical Device directive 93/42/EEC but, being a developer's version, it's not certified.

## Content

<b>BioLib library</b> .....	4
<b>Class:</b> data parser.....	5
<b>Methods</b> .....	5
<b>Events</b> .....	7
<b>Class:</b> accelerometer data.....	9
<b>Class:</b> Qrs data.....	10
<b>Class:</b> utils.....	11
<b>Functions</b> .....	11
<b>References</b> .....	13
<b>Control versions</b> .....	14



## BioLib library



public class **Data**  
Member of [BioLib](#)

public class **DataACC**  
Member of [BioLib.Data](#)

**Summary:**  
Accelerometer data.

public class **Qrs**  
Member of [BioLib.Data](#)

**Summary:**  
Qrs data.

public class **Utils**  
Member of [BioLib](#)

**Summary:**  
Utils functions.

## Class: data parser

- Data(string, BioLib.Data.LEADTOANALYSE, BioLib.Data.NUMBERPACKETPERSECOND)
- GetEcgOffset()
- GetNumberOfLeads()
- GetSampleFrequency()
- GetVersion()
- SetData(byte[], int)
- ⚡ accConfEventHandler
- ⚡ accEventHandler
- ⚡ batteryEventHandler
- ⚡ deviceIdEventHandler
- ⚡ ecgEventHandler
- ⚡ firmwareVersionEventHandler
- ⚡ pushbuttonEventHandler
- ⚡ qrsEventHandler
- ⚡ radioEventHandler
- ⚡ rtcEventHandler
- ⚡ sdcardEventHandler
- ⚡ timestampEventHandler

## Methods

public **Data**(string *applicationPath*, [BioLib.Data.LEADTOANALYSE](#) *lead*, [BioLib.Data.NUMBERPACKETPERSECOND](#) *nPacket*)

Member of [BioLib.Data](#)

### Summary:

Constructor.

### Parameters:

*applicationPath*: Application path

*lead*: lead use to detect QRS

*nPacket*: number of packet ECG per second (1 – 500 samples ECG / sec., 5 – 100 samples ECG / sec. or 10 – 50 samples ECG / sec.)

public [System.Collections.Generic.List<string>](#) **GetNameOfLeads**()

Member of [BioLib.Data](#)

### Summary:

Get name of leads.

### Returns:

Name(s) of lead(s)

public [int](#) **GetNumberOfLeads()**  
Member of [BioLib.Data](#)

**Summary:**

Get number of leads.

**Returns:**

Number of leads

public [int](#) **GetSampleFrequency()**  
Member of [BioLib.Data](#)

**Summary:**

Get ECG sample frequency.

**Returns:**

ECG sample frequency (Hz)

public [bool](#) **SetData([byte\[\]](#) data, [int](#) nBytes)**  
Member of [BioLib.Data](#)

**Summary:**

Set data to parser.

**Parameters:**

*data*: Data to parser

*nBytes*: Number of bytes in buffer

**Returns:**

True, if no errors occurred

public [string](#) **GetVersion()**  
Member of [BioLib.Data](#)

**Summary:**

Get BioLib library version.

**Returns:**

Version of library.

## Events

public event [BioLib.Data.AccEventHandler](#) **accEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.BatteryLevelEventHandler](#) **batteryEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.DeviceIdEventHandler](#) **deviceIdEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.EcgEventHandler](#) **ecgEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.PushbuttonEventHandler](#) **pushbuttonEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.QrsDetectedEventHandler](#) **qrsEventHandler**  
Member of [BioLib.Data](#)

QRS detector is based on the algorithm of Pan and Tompkins [1] and was used MIT-BIH database for validate results [2].

public event [BioLib.Data.RadioEventHandler](#) **radioEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.RTCEventHandler](#) **rtcEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.SdCardStatusEventHandler](#) **sdcardEventHandler**  
Member of [BioLib.Data](#)


public event [BioLib.Data.TimestampEventHandler](#) **timestampEventHandler**  
Member of [BioLib.Data](#)

public event [BioLib.Data.FirmwareVersionEventHandler](#) **firmwareVersionEventHandler**  
Member of [BioLib.Data](#)





## Class: accelerometer data



public class **DataACC**  
Member of [BioLib.Data](#)

### Summary:

Tri-axial accelerometer data (X,Y,Z).

public [int](#) **Position** { set; get; }  
Member of [BioLib.Data.DataACC](#)

### Summary:

Position in ecg samples offset.

public [sbyte](#) **X** { set; get; }  
Member of [BioLib.Data.DataACC](#)

### Summary:

Axis xx.

public [sbyte](#) **Y** { set; get; }  
Member of [BioLib.Data.DataACC](#)

### Summary:


Axis yy.

public [sbyte](#) **Z** { set; get; }  
Member of [BioLib.Data.DataACC](#)

### Summary:

Axis zz.

## Class: Qrs data



public class **Qrs**  
Member of [BioLib.Data](#)

**Summary:**  
Qrs data.

public [short](#) **bpmi** { set; get; }  
Member of [BioLib.Data.Qrs](#)

**Summary:**  
Bpm instantaneous (bpm).

public [int](#) **position** { set; get; }  
Member of [BioLib.Data.Qrs](#)

**Summary:**  
Qrs position (samples offset).

public [int](#) **rr** { set; get; }  
Member of [BioLib.Data.Qrs](#)

**Summary:**  
Beat-by-beat R-R (ms).

## Class: utils

- ⇒ BCDtoDEC(int)
- ⇒ DECtoBCD(int)
- ⇒ Get\_CRC(byte[], int)
- ⇒ GetPulse(float)
- ⇒ GetRR(int, int)

## Functions

public static [int](#) **BCDtoDEC**([int](#) bcdValue)  
Member of [BioLib.Utils](#)

### Summary:

Convert BCD format in Decimal format.

### Parameters:

*bcdValue*: BCD value

public static [byte](#) **DECtoBCD**([int](#) value)  
Member of [BioLib.Utils](#)

### Summary:

Convert decimal format in BCD format.

### Parameters:

*value*: decimal value

public static [byte](#) **Get\_CRC**([byte\[\]](#) data, [int](#) nBytes)  
Member of [BioLib.Utils](#)

### Summary:

Calculate CRC of stream data (send to device).

### Parameters:

*data*: stream data

*nBytes*: number of bytes of stream

### Returns:

value CRC

public static [float](#) **GetPulse**([float](#) rr)  
Member of [BioLib.Utils](#)

**Summary:**

Get BPM instantaneous (bpm).

**Parameters:**

*rr*: R-R (ms)

public static [float](#) **GetRR**([int](#) rr, [int](#) SampleFrequency)  
Member of [BioLib.Utils](#)


**Summary:**

Get R-R (ms).

**Parameters:**

*rr*: R-R (samples)

*SampleFrequency*: Sample Frequency (Hz)



## References

- [1] Pan J and Tompkins WJ. A Real-Time QRS Detection Algorithm. IEEE Transactions on Biomedical Engineering 32(3):230-236, 1985
- [2] MIT-BIH Arrhythmia Database: <http://www.physionet.org/physiobank/database/mitdb/>



## Control versions

Version	Date	Change log
1.0.02	30-04-2013	<i>Get device ID</i>
		<i>Send radio event to device</i>
1.0.03	19-07-2013	<i>New method to send radio-event to device</i>
		<i>A new too l(InfoExporter.exe) for export data to Excel and Matlab</i>
1.0.07	18-03-2015	<i>Get firmware version from device VitalJacket.</i>

